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**Directions of use**

# **C-Cement Flow**

**chemisch & lichthärtender Komposite-Zement**

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## 1. Specification of medical device product

C-Cement Flow

C-Cement Flow is a dual cure flowable Microhybrid-Composite for fixation of veneers, inlays, crowns and bridges. C-Cement Flow based on base and catalyst in an auto mix syringe. C-Cement Flow is X-ray opaque and has high filler content. This guaranties high colour stability and a high abrasion resistance.

## 2. Intrigents

*filler:*

The filler of C-Cement Flow is based on SiO<sub>2</sub> and other fine glasses (Barium-Aluminium-Bor-Silikatglass).

*Acrylic matrix:*

The organice matrix of C-Cement Flow consist of:

2,2-Bis-4(2hydroxy-3-methacryloxy-propyloxy)-phenyl-propan (Bis-GMA)

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3,6-Dioxaoctamethylendimethacrylat (TEGDMA).

## 3. Indicationen

For fixation of Inlays, Onlays, veneers and crown & bridge works

## 4. side attention

not appliance.

## 5. Working instruction

C-Cement Flow is dual curing flow able resin cement in a atomix-syringe.

The automix-syringe is ready for direct application after fixation of the mixing tip on the top.

The application time (at 23°C) is 2 min after auto mixing. The chemical polymerisation starts after mixing process within 5 minutes (2 min. application + 3 min polymerization). C-Cement can be cured with dental light, e.a. Megalux CS in 40 sec. Without light (e.a. under crowns, etc.) the material will be chemical cured in 3 minutes after application.

The cavity can be prepared like the normal adhesiv technique. The under filling must be resisted against acid. The prepared oversize of the teeth must be cleaned with cleaning instruments.

The cavity must be acid with C-Cid etching gel for 30 sec. After this time, C-Cid must be cleaned with water. Please dry the cavity with oil-free air. Apply C-Bond in a tine layer on the whole cavity. After application the cavity will be ready for C-Cement Flow.

## 6. Shelf life

2 years

## 7. Storage advice

cool and dry storage (+5°C till +8°C) ,  
please not over +25°C.

## 8. Time of information

January 2011

## 9. Name and address of manufacture

MEGADENTA Dentalprodukte GmbH  
Carl-Eschebach-Strasse 1 A  
D-01454 Radeberg  
Germany