## Directions of use

### 1. Information

Megafill Flow, flowable light-curing Microhybridcomposite  
Shades: A1 / A2 / A3 / A3,5 / B2

### 2. Composition

The inorganic filler particles comprise barium aluminium bore glass and highly dispersed silicon dioxide.

The monomer matrix contains:

2,2-Bis-4(2-hydroxy-3-methacryloxy-propyloxy)-phenyl-propan;

7,9,9-trimethyl-4,13-dioxo-3,14-dioxo-5,12-diaza-hexadecan-

1,16-di-yl-bis-methacrylate,

1,4-Butandiol dimethacrylate.

### 3. Indications

Class V restorations (cervical caries, root erosion, wedge-shaped defects),

small posterior restorations (classes I and II),

anterior restorations (classes III and IV),

restorative therapy for mini-cavities,

preventive resin restorations in premolars and molars,

adhesive cementation of ceramic and composite restorations.

### 4. Side Effects

Placement of Megafill Flow is contraindicated if the patient is known to be allergic to any of the ingredients in Megafill Flow.

Irritations resulting from direct contact with dentin in proximity to pulp cannot be ruled out. In such cases the dentin should be shielded with a calcium hydroxide compound and a subfilling (e. g. glass ionomer cement)

### 5. Interaction

Application of eugenol containing material, specifically base liners, is not indicated, it will impair the polymerization.

### 6. Warning

Since Megafill Flow cures with light, it is sensitive to blue light (operating light, sun light). For that reason bright light during application of Megafill Flow should be avoided.

Megafill Flow should not come in contact with skin / mucous membrane or eyes.

## 7. Application

#### - Shade Selection

The tooth and the adjacent distal surfaces has to be cleaned with a bristle brush and prophylactic paste. Shade matching should preferably be carried out in daylight using the MEGADENTA Shade Guide. The shades are corresponding to the Vita-Lumin-Vacuum® Shade Guide.

#### - Cavity Preparation

If the working field cannot be established dry the application will be contraindicated. The cavity is prepared and dried (rubber dam) in accordance with the general rules of composite filling treatment and unprepared dentine is cleaned with a cleaning paste. Caries-free cervical defects are not prepared, but merely cleaned with pumice and a suitable cleaning paste using a rubber cup or rotary brush. The prepared cavity is now rinsed with water and dried with oil free air. Deeper lying areas of the cavity can be filled with a suitable under filling material.

#### - Base Liner

To protect exposed dentine the application of a glass ionomer liner, phosphate cement or calcium hydroxide base liner (deep cavities) is recommended, thus is reducing the danger of secondary caries at the margins.

#### - Acid Etching

The cleaned enamel edges are etched with C-Cid. C-Cid is applied with the supplied brush. After 30 – 45 s C-Cid is applied on the dentine allowed to take effect for another 10 – 15 s and then thoroughly rinsed with water. Excessive drying should be avoided in this case, as this would cause the exposed collagen fibre to collapse.

#### - Application of the Primer

The Primer has now to be used after the recommendations of the producer.

Do not rinse off the dentin primer!

The cavity is dried thoroughly with oil and water free air or until the dentine surface appears dull and dry. It is essential that the primed dentine and the etched enamel are dry and contaminant free for the bond application.

C-Bond is applied onto the primed dentine in the dried cavity and the etched enamel, blown out thinly with oil free air and then cured for 10 s with a blue light source (e.g. Megalux CS). Following the polymerization of C-Bond Megafill Flow can be applied.

#### - Application of Megafill Flow

The material is applied directly into the cavity. To avoid an oxygen inhibited layer, a coat of oxygen-tight glycerine gel can now be applied along the restoration margins. Make sure that the glycerine gel is not mixed with the Megafill Flow.

After completed contouring the restoration is cured for 40 s with Megalux CS unit or an other halogen light unit. The beam has to be directed to all surfaces of the restoration.

Restorations deeper than 2,5 mm should be filled in layers, each layer have to be cured separately (multi-phase technique).

#### - Shaping and Finishing

The final forming and the removal of excess composite are carried out with rotating finishing instruments. For contouring and finishing fine and ultra fine finishing diamonds are used. For polishing flexible discs and polishing strips are suitable. All shaping and finishing have to be carried out under tooth preserving and with water-cooling. A sealing of the filling with C-Bond is recommendable.

### 8. Shelf life

3 Years

### 9. Precautions and Storage

Not to be stored above 25 °C.

Storage in refrigerator (+5 to +8 °C) prolongs shelf life.

### 10. Deliver units

Megafill Flow single syringe 2,0 g

Megafill Flow assortment with 4 syringes + accessories

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MEGADENTA Dentalprodukte GmbH
Wir, MEGADENTA Dentalprodukte GmbH, 01454 Radeberg, erklären in alleiniger Verantwortung, dass das Produkt
Megadenta Dentalprodukte GmbH