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Megaseal Universal

properties and advantages:

- **easy and fast application**
- **high fitting accuracy**
- **high gloss oversize**
- **cureable in all dental light units**

1. Product description

Megaseal Universal is a light curing lacquer for the coating of resin parts of acrylic dentures as well for temporaries in C+B technique.

An abrasion-resistant and polish like surface can be produced.

2. Application

Megaseal Universal can be applied with a brush directly onto the oversize of the denture or temporary. For safe adhesion the surface should be free of grease and dry.

The layer should be thin. If a thicker layer is required, the application should be in 2 steps with each separate curing process.

The hardening / curing process should be effected in a light-polymerization unit (e.a. Megalight) with a wave length of 320-500 nm. Based on the intensity of the unit, the curing time should be 3-5 minutes.

3. Important advice

Attention! Megaseal Universal contains Methylmethacrylate ! It is flammable liquid. Keep away from fire or heat. Do not inhale and do not get in contact with skin or eyes.

Contains: Methylmethacrylate, Polyacrylate, Photoinitiator, stabilizer