Gebrauchsanweisung
Instruction for use
Mode d'emploi
Istruzioni d'uso
Mode de empleo

Megatray

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Instruction for use

(GB/IR/US/CA)

1. Name and description of the medical device

MEGATRAY consists of a matrix of polyfunctional methacrylates and inorganic fillers. The resin is slightly transparent, blue or pink. The thickness of the MEGATRAY tray material is 2.2 mm. The thickness of MEGATRAY base plate is 1.4 mm.

MEGATRAY has

- less shrinkage,
- resistance to water absorption
- high dimensional stability.

2. Indication

Light-cure impression tray material for the fabrication of custom trays, bite plates and base plates for needle-point tracing.

3. Application

The plaster model must be insulated with alginate insulation (MEGA-I) before using MEGATRAY. The cover-vlies from the wafer have to be removed, than MEGATRAY is moulded onto the model. The tray can be cut to the desired rim. Use excess material to form the handle. The wafer is polymerized in a UVA- or blue-light unit (e.a. Megalight ST, Megalight Mini) for 5 minutes. The wafer is taken off the model and the reverse side of the tray is cured for another 5 minutes. After polymerising the inhibition layer caused by the airoxygen is removed by a tissue moistened with an alcoholic cleaner (e.a. Megaclean) easily. Use a dust extractor for trimming with carbide burs.

4. Indication of risk

MEGATRAY contains monomers for polymerization (e.g. methacrylic esters) which can cause skin allergies, especially with sensitive patients. Do not use this product where allergic reactions are observed, or for patients with methacrylate/acrylate-based resin allergy history.

5. Shelflife

24 month

6. Storage

Not to be stored above 25°C / 77°F

7. Package

Standardpackage:

100 wafers upper slightly transparent	Art.-no. 1 MYEVOK-TR
100 wafers upper pink	Art.-no. 1 MYEVOK-RS
100 wafers upper blue	Art.-no. 1 MYEVOK-BL
50 wafers upper slightly transparent	Art.-no. 1 MYDVOK-TR
50 wafers upper pink	Art.-no. 1 MYDVOK-RS
50 wafers upper blue	Art.-no. 1 MYDVOK-BL

Accessories:

Megaclean	alcoholic cleaning agent
Megatray Adhesiv	conector to silicon mat.
Mega-I	alginat-insolation
Megaseal UNIVERSAL	sealed agent, tissue glaze

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MEGADENTA Dentalprodukte GmbH