Megafill Bio

[de]	Gebrauchsanweisung
[bg]	Инструкции за употреба
[cs]	Návod k použití
[en]	Instructions for use
[es]	Instrucciones de uso
[fr]	Mode d'emploi
[hr]	Upute za uporabu
[hu]	Használati utasítás
[it]	Istruzioni per l'uso
[lt]	Naudojimo instrukcija
[lv]	Lietošanas pamācība
[pl]	Instrukcja używania
[ro]	Instrucțiuni de utilizare
[ru]	Инструкция по применению
[sk]	Návod na použitie
[sv]	Bruksanvisning

MEGADENTA
Dentalprodukte

MEGADENTA Dentalprodukte GmbH
Carl-Eschebach-Str. 1A, D-01454 Radeberg, Germany
Telefon +49(0)3528 453-0, Fax +49(0)3528 453-21
E-Mail: info@megadenta.de www.megadenta.de

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Read and keep before using the product! Medical device for professionals only

1 General safety instructions

Avoid contact of uncured material with skin/mucous membrane and eyes. It may be slightly irritating and in rare cases may cause sensitization to methacrylate.

Use protective equipment: Safety glasses, suitable protective gloves with protection function to Methacrylates (e.g. Polyethylene-gloves, 4h-gloves).

After skin contact: Wash immediately with plenty of soap and water. Medical treatment necessary if symptoms occur (irritation of the skin).

After eye contact: Immediately flush with running water for 10 to 15 minutes with the eyelid open and consult an ophthalmologist.

If swallowed: Rinse your mouth immediately and drink plenty of water.

In the course of a filling's service life, which may be several years, the possibility of chemicals/substances dissolving out of the resin matrix cannot be completely excluded. This depends on various factors (e.g. nutritional habits of the patient). The patient must be informed about the residual risk of genotoxic, cytotoxic and/or sub-acute systemic toxic effects, with a special focus on pregnant women, nursing mothers and children.

The nanomaterial contained in the polymerized state does not cause any increased risk regarded to cell-permeable nanoparticles, since the polymerization results in a solid compound.

The publicly available Summary Safety and Clinical Performance Report (SSCP) according to Article 32 of Regulation (EU) 2017/745 (MDR) can be found in the EUDAMED database at the following link: <https://ec.europa.eu/tools/eudamed/#/screen/search-device>

All serious incidents related to the product must be reported to the manufacturer and the competent authority of the member state where the user and/or patient is established. Please direct any feedback regarding the product to: ***reclamation@megadenta.de***

2 Description

2.1 General description

Megafill Bio is a polymer-based restorative material class 2, group 1 and 2, classified according to DIN EN ISO 4049:2019.

The material is a paste-like, modellable, light-curing, low-allergen universal microhybrid composite for anterior and posterior restorations in the shades*: Enamel: A1, A2, A3, A3.5, B2, B3, C2, D3; Dentin: A2, A3 (*according to Vita Classic shade system of Vita Zahnfabrik, Bad Säckingen, Germany)

Megafill Bio is suitable for the restoration of occlusal surfaces.

The dental fillings made of Megafill Bio are radiopaque (2,2 Al equivalent).

Note: Aluminum has a radiopacity equivalent to dentin. Therefore, 1 mm of a material that has a radiopacity equivalent to 1 mm of aluminum has a radiopacity equivalent to 1 mm of dentin and 2 mm of aluminum is equivalent to enamel.

2.2 Composition

Silanized glass, UDMA, pyrogenic silica, silanized SiO₂, BDDMA, initiators/stabilizers/pigments each < 0.1 %; mass fraction of monomers ~ 23 %, volume fraction of fillers ~ 59 %, particle size range of inorganic fillers 0.2-0.7 µm.

3 Intended use

3.1 Warnings

Megafill Bio is light-sensitive, so that polymerization may already be triggered by exposure to daylight or by ambient light from the dental treatment lamp. The material should therefore be removed from the primary packaging only immediately before application and the packaging resealed directly. In case of damage to the primary packaging (syringe), the product must be disposed according to 3.12.

After application, the material in the syringe should be relaxed by turning back the plunger.

Single-dose packaging (Minifills) should be disposed of after completion of the application (see 3.12).

For hygienic reasons (see 4.6), they are intended for single use only.

3.2 Intended purpose / intended use

Megafill Bio is intended for dental use for the restoration of dental defects:

- filling of dental cavities of classes I, II, III, IV, V
- composite inlays
- fixation and repair of composite and ceramic restorations
- direct composite veneers
- splinting of loosened teeth
- Undercut block out

The user is obliged to use the material in accordance with the intended purpose specified by the manufacturer MEGADENTA Dentalprodukte GmbH. The manufacturer accepts no liability for damage resulting from use for other purposes.

3.3 Clinical benefit

The clinical benefit of Megafill Bio is the improvement or restoration of functionality and/or esthetics of teeth.

3.4 Indications

- Class I, II, III, IV, V dental cavities
- defects on composite or ceramic restorations
- Discoloration of anterior teeth
- Loosened teeth
- Undercuts

3.5 Contraindication

- insufficient moisture control of the working area
- known hypersensitivity to the substances listed in the composition (especially methacrylate)
- eugenol-containing underfillings (see 3.9)

3.6 User profile

Megafill Bio is a product for professional use by dentists. Application is intended only in a dental practice. Special training is not required.

3.7 Patient group

The use of Megafill Bio is not limited to any particular patient group. Contraindications and safety instructions must be observed.

3.8 Undesirable side effects

Megafill Bio contains methacrylate. These may cause allergic skin reactions. In rare cases, sensitization is possible. Pulpal irritation cannot be excluded in the case of cavities close to the pulp. A protective base/ underfilling is recommended in such cases.

3.9 Undesirable interactions

Megafill Bio must not be used together with eugenol-containing underfillings, as phenolic substances inhibit polymerization.

3.10 Special storage instructions

Store protected from light and at 4-25 °C.

Storage outside the specified conditions may result in damage to the product, e.g. segregation, worse mechanical performance and reduced applicability, premature curing. It can then no longer be applied as intended and must be disposed of in accordance with 3.12.

3.11 Shelf life

3 years. Do not use after the expiry date stated on the label!

3.12 Waste treatment

Do not allow product or primary packaging which has not been emptied of all residues to reach sewage system, water, subsoil or soil. Disposal according to official regulations.

Non-contaminated and completely empty packaging as well as fully cured products can be disposed of with normal municipal waste (household waste). Observe waste separation.

4 Application

4.1 Shade determination

Clean the teeth with oil-free material before determining the shade. The color is determined on the still moistened tooth using the Vita Classic color system.

4.2 Drying

Thoroughly dry the surfaces to be treated. The use of dental dams is recommended.

4.3 Cavity preparation

The cavity preparation is carried out according to the rules of the adhesive filling technique. In the anterior region, bevel the enamel margins; in the posterior region, only slightly break or round off the sharp enamel edges, always sparingly and conserving the substance. Then remove all residues and dry the cavity.

4.4 Underfilling

If an underfilling is to be placed (e.g. for cavities close to the pulp), make sure that the material is acid-resistant and eugenol-free.

4.5 Adhesive base preparation

Prepare the bonding surfaces and apply the bonding agent according to the etching technique and in accordance with the instructions for use of the products being used.

Megafill Bio is compatible with all fixing materials from MEGADENTA Dentalprodukte GmbH.

4.6 Application

- Allow material to acclimatize to room temperature before application to minimize cold irritation.
- Apply the filling material into the prepared cavity using a metal or plastic spatula or a ball tamper. When using Minifills, the product can be applied directly from the packing into the cavity. A Minifill is only intended for the treatment of a single patient and must be disposed of after completion of the treatment according to item 3.12 to avoid cross-contamination.
- The layer thickness must not exceed 2 mm. In the case of voluminous fillings, the multi-layer technique must be used and each layer must be cured individually.
- Curing is carried out with a commercially available blue light source, e.g. Megalux CS (emission wavelength range 400-500 nm, irradiance 600-800 mW/cm²).
- The exposure time with the recommended device is uniformly 40 s for all colors. The polymerization depth after this irradiation time is 2 mm.
- The instructions for use of the device used must be observed!

4.7 Finishing/control/polishing

Caution! Avoid inhaling the dust during grinding/polishing, use mouth and nose protection.

After the prescribed exposure time, the filling can be processed immediately with the usual rotary instruments (carbide drill, diamond, finisher, polisher). Check occlusion and articulation and grind in if necessary. Polishing the filling improves the marginal seal and the esthetics of the filling.

5 Delivery units

Single syringe 4.5 g; Assortment packs; Minifills ( 20 x Minifill á 250 mg)

6 Symbols used in product labeling

	Manufacturer
	Medical device
	Item number
	Batch
	Usable until
	Protect from sunlight
	Temperature limit
	Do not reuse
	Contains nano materials
	Observe instructions for use (The

	round ISO 7010-M002 symbol on the packaging is blue)
	European Conformity Mark

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Megafill Bio är kompatibelt med alla fästmaterial från MEGADENTA Dentalprodukte GmbH.