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## **Instructions for use**

### **1. Product description**

Bloxx cast is a light curing resin for the generative production of burn-out objects for casting process, application for DLP printers with a 385 nm UV light source.

### **2. Composition**

Combination of acrylic resin (Methacrylate base) with Photoinitiators and stabilisators

### **3. Indication**

Generative production of burn-out objects for casting process in dental and jewelry

### **4. Contraindications**

Skin contact should be avoided in uncured acrylic resin.

### **5. Side effects**

In individual cases, contact allergy may occur.

### **6. Interactions**

Unknown.

### **7. Hazard Information**

Read Safety data sheets first!

Avoid skin contact when working with methyl methacrylates, direct skin and eye contact. (Use personal protective equipment: protective gloves / eye protection). If swallowed, seek medical advice immediately. Ensure adequate ventilation during processing. Do not inhale vapors.

## **8. Handling**

### ***Preparations***

Storage between 15°C-28°C, do not shake before using

### ***Processing***

Bloxx cast has been optimised for the DLP printing system. Following the printing process, the printed objects must be cleaned, dried and post-exposed with a light source, e.a. Otofash G171 xenon photoflash unit to guarantee the requisite product characteristics.

- layer thickness: 50 µm
- ensure the material is free of bubbles before printing process
- dripping time of approx. 10 minutes once printing recommended
- finish the printed objects without delay so as to avoid distortion of their shape and guarantee the required precision

### **Post-processing:**

#### *- Cleaning*

We recommend a solvent-resistant cleaning mat for cleaning. Detach the printed objects carefully from the building platform. Cleaning should be performed with Isopropanol ( $\geq 98\%$  purity) in an unheated ultrasound bath. The operating instructions from the manufacturers must be observed and followed. The printed objects should be cleaned in two stages.

(1) Pre-cleaning: Pre-clean the printed object for 2 minutes in a reusable Isopropanol ultrasound bath. Note: The bath's cleaning efficacy decreases the more it is used. If the cleaning performance deteriorates, the bath should be replaced.

(2) Main cleaning process: Then clean the printed objects for a further 2 minutes in a fresh Isopropanol ultrasound bath. Should there still be resin residue on the printed object after the main cleaning process, the printed object can be cleaned again for max. 1 minute in the Isopropanol.

#### *- Finishing*

Remove and sand the support structures carefully using a rotary instrument and without exerting pressure prior to the post-exposure. If no more resin residues are visible on the printed object, the printed object should be dried carefully using compressed air.

#### *- Post-exposure*

The post-exposure is performed twice with an Otofash G171 xenon photoflash unit with 2,000 flashes each site or a similar light source, protective gas atmosphere is recommended. Following the first 2,000 flashes, allow a cooling-off period of at least 2 minutes with the lid open. Insufficient cooling can distort the shape of the printed object. Once the cooling-off period has elapsed, the printed objects should be turned and light-cured with a further 2,000 flashes.

Embedding and casting: Please use phosphate-bonded bedding compounds for production of casting mould only. Follow the instruction for use of the bedding compound and device manufactures.

## **9. Shelf life**

2 years

## **10. Precaution and storage**

Storage between 15°C and 28 °C

## **11. Delivery units**

Liquid: 1000 g shade: red

## **12. Waste treatment**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation. Handle contaminated packages in the same way as the substance itself.

**13. Date information prepared**

*2020-01-31*

**14. Manufacturer address**

MEGADENTA Dentalprodukte GmbH  
Carl-Eschebach-Str. 1A / D-01454 Radeberg, Germany

---

*1100 0120 541*