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Instructions for use

1. Product description

Bloxx model is a light curing acrylic resin for the generative production of dental models for DLP printers with a 385 nm UV light source.

2. Composition

Combination of acrylic resin (Methacrylate base) with Photoinitiators, stabilisators and fillers

3. Indication

Generative production of dental working and presentation models

4. Contraindications

Skin contact should be avoided in uncured acrylic resin.

5. Side effects

In individual cases, contact allergy may occur.

6. Interactions

Unknown.

7. Hazard Information

Read Safety data sheets first!

Avoid skin contact when working with methyl methacrylates, direct skin and eye contact. (Use personal protective equipment: protective gloves / eye protection). If swallowed, seek medical advice immediately. Ensure adequate ventilation during processing. Do not inhale vapors.

8. Handling

Preparations

Storage between 15°C-28°C, shake well before using

Processing

Bloxx model has been optimised for the DLP printing system. Following the printing process, the printed objects must be cleaned, dried and post-exposed with a light source, e.a. Otofash G171 xenon photoflash unit to guarantee the requisite product characteristics.

- layer thickness: 50 µm
- ensure the material is free of bubbles before printing process
- recommended wall thickness if the model is not filled: 3 mm
- for precision work, place models horizontally and cores vertically on the building platform
- for more efficient printing, deep drawn and antagonist models can be positioned vertically on the building platform.
- dripping time of approx. 10 minutes once printing recommended
- finish the printed objects without delay so as to avoid distortion of their shape and guarantee the required precision

Post-processing:

- Cleaning

We recommend a solvent-resistant cleaning mat for cleaning. Detach the printed objects carefully from the building platform. Cleaning should be performed with Isopropanol ($\geq 98\%$ purity) in an unheated ultrasound bath. The operating instructions from the manufacturers must be observed and followed. The printed objects should be cleaned in two stages.

(1) Pre-cleaning: Pre-clean the printed object for 2 minutes in a reusable Isopropanol ultrasound bath. Note: The bath's cleaning efficacy decreases the more it is used. If the cleaning performance deteriorates, the bath should be replaced.

(2) Main cleaning process: Then clean the printed objects for a further 2 minutes in a fresh Isopropanol ultrasound bath. Should there still be resin residue on the printed object after the main cleaning process, the printed object can be cleaned again for max. 1 minute in the Isopropanol.

- Finishing

Remove and sand the support structures carefully using a rotary instrument and without exerting pressure prior to the post-exposure. If no more resin residues are visible on the printed object, the printed object should be dried carefully using compressed air.

- Post-exposure

The post-exposure is performed with an Otofash G171 xenon photoflash unit with 2 x 2,000 flashes or a similar light source, protective gas atmosphere is recommended. Following the first 2,000 flashes, allow a cooling-off period of at least 2 minutes with the lid open. Insufficient cooling can distort the shape of the printed object. Once the cooling-off period has elapsed, the printed objects should be turned and light-cured with a further 2,000 flashes.

9. Shelf life

2 years

10. Precaution and storage

Storage between 15°C and 28 °C

11. Delivery units

Liquid: 1000 g shade: beige

12. Waste treatment

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation. Handle contaminated packages in the same way as the substance itself.

13. Date information prepared

2020-01-31

14. Manufacturer and address

MEGADENTA Dentalprodukte GmbH
Carl-Eschebach-Str. 1A / D-01454 Radeberg, Germany

1100 0120 540