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
Instructions for use

1. Information

C-Seal is a light-curing pit and fissure sealant.

Variants: transparent or white paste with low viscosity

2. Composition

C-Seal contains highly dispersed silicon dioxide and Silicate glass as inorganic fillers (47.7 wt.-%) and Fluorine as Aminfluoride (0.15 wt.-%). The monomers are Bis-GMA, Bis-EMA and Butandioldimethacrylate (49.8 wt.-%).

3. Indication

Usage in dentistry for prophylactic treatment for sealing fissures, for caries prophylaxis, for restoring small carious defects and for sealing fillings made of plastic or cement (moisture protection)

4. Contraindication

- If a dry working field cannot be established.
- If the patient is known to be allergic to any of the ingredients.

5. Side effects

The unpolymerized material may cause contact allergy.

6. Other interactions

C-Seal may not be used together with underfillings containing eugenol because phenolic substances inhibit the polymerisation process.

7. Warning

C-Seal is light sensitive as polymerisation can happen with day-light. Therefore the light exposure should be limited (take material from syringe just before application, close syringe directly after that).

8. Hazard Information

Read Safety data sheet first!

General information

Polymerized material is non-hazardous.

After inhalation

Provide fresh air. When in doubt or if symptoms are observed, get medical advice.

After contact with skin

After contact with skin, wash immediately with polyethylene glycol, followed by plenty of water. Take off immediately all contaminated clothing and wash it before reuse. Medical treatment necessary.

After contact with eyes

In case of contact with eyes flush immediately with plenty of flowing water for 10 to 15 minutes holding eyelids apart and consult an ophthalmologist.

After ingestion

Rinse mouth immediately and drink plenty of water.

Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

9. Application

-Cleaning and drying the teeth:

- Before etching, the tooth surface has to be thoroughly cleaned with an oil free polishing paste or a pumice powder.
- Dry the surfaces.
- Use of kofferdam is recommended.

-Enamel etching:

- Use the supplied brush to etch the cleaned enamel with C-Cid.
- Rinse C-Cid with water after 20-40s thoroughly and dry carefully with oil- and water-free air. The etched and dried enamel must show a smooth, chalky surface.
- The etched enamel must not be contaminated by saliva or blood since this would lead to a deposit of protein, which would inactivate the etched enamel surface. If this happens, repeat the etching process.

-Application:

- Apply C-Seal directly with the needle into the fissure. Direction of application is recommended from the margin into the fissure.
- Wait for 15-20s for penetration of sealant.
- Remove any air bubbles with the needle or with a suitable instrument.
- Polymerize with a light source (e.g. Megalux) for 20 – 40 s.
- Clean the sealed surface.
- In case of “white” C-Seal the sealed parts of the tooth are easy detectable by its white shade.

10. Shelf life

2 years

11. Precaution and Storage

We recommend storage in a fridge (+5 °C [41 °F] to +8 °C [46 °F]); do not exceed 25 °C [77 °F].

12. Delivery units

single syringe C-Seal á 1.8 g white: REF 2CSEAL-OP transparent: REF 2CSEAL-TR
assortment 4 x C-Seal á 1.8 g REF 2CSS
Accessories: brush holder, brush inserts, Luer-Lock cannulas

13. Waste treatment

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

Keep out of the reach of children!
For use in dentistry only!

14. Date information prepared

2019-04-05, Rev. 4

15. Manufacturer and Sales

MEGADENTA Dentalprodukte GmbH
Carl-Eschebach-Strasse 1 A
D-01454 Radeberg, Germany

1340 0419 529