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### 1. Information

Mega-A

Cold curing resin for repairs and relining of total and partial dentures

acc. EN ISO 20795-1, class 1, type 2

Shades: pink, blue-pink veined and clear

### 2. Composition

Powder/liquid cold cure denture resin

Main component powder: polymethyl methacrylate

Component liquid: methyl methacrylate, dimethacrylate

### 3. Areas of application

Resin for functional reconstruction, repairing, extending, supplementing as well as partial and total re-lining of dentures

### 4. Contraindications

Skin contact should be avoided in patients with known hypersensitivity to methacrylates.

### 5. Side effects

In individual cases, contact allergy may occur.

### 6. Interactions

Unknown.

### 7. Hazard Information

Read Safety data sheets first!

Mega-A liquid is highly flammable, the vapors may form explosive mixtures with air. Remove all sources of ignition. Close bottle immediately after usage.

Avoid skin contact when working with methyl methacrylates, direct skin and eye contact. (Use personal protective equipment: protective gloves / eye protection). If swallowed, seek medical advice immediately. Ensure adequate ventilation during processing. Do not inhale vapors.

### **Mega-A Powder:**

#### *AFTER INHALATION*

*Provide fresh air. If breathing is irregular or stopped, administer artificial respiration. Medical treatment necessary.*

#### *AFTER CONTACT WITH SKIN*

*After contact with skin, wash immediately with polyethylene glycol, followed by plenty of water. Take off immediately all contaminated clothing and wash it before reuse. If skin irritation occurs: Get medical advice/attention.*

*AFTER CONTACT WITH EYES*

*Rinse immediately carefully and thoroughly with eye-bath or water.*

*AFTER INGESTION*

*Rinse mouth immediately and drink plenty of water.*

***Mega-A Liquid:***

*AFTER CONTACT WITH SKIN*

*After contact with skin, wash immediately with polyethylene glycol, followed by plenty of water. Take off immediately all contaminated clothing and wash it before reuse. Medical treatment necessary.*

*AFTER CONTACT WITH EYES*

*After contact with the eyes, rinse with water with the eyelids open for a sufficient length of time, then consult an ophthalmologist immediately.*

*AFTER INGESTION*

*Rinse mouth immediately and drink plenty of water.*

## **8. Handling**

### *- Preparations*

For embedding the prosthesis set up in wax, a semi-hard plaster is sufficient. After the wax has been scalded, the plaster is coated with MEGA-I or another alginate insulating agent.

### *- Bond to acrylic- and ceramic teeth*

Roughen the basal surface of the acrylic teeth before bonding. Any contamination by wax and isolating agent must be avoided. Mechanical retention is recommended. With ceramic teeth, sufficient mechanical retention and, if necessary, additional treatment (silanization) must be ensured. Please observe the processing instructions of the tooth manufacturer.

### *- Dosage and Handling*

Powder and liquid are intensively mixed in a ratio of **5 g powder : 3 ml liquid** within 20 s. The liquid is placed in the mixing vessel and the powder is interspersed until a low liquid level remains. The droppable liquid dough can be processed immediately. Apply MEGA-A evenly. After reaching the plastic consistency, do not touch and model again.

### *- Polymerization*

The Curing time at 21 °C is approximately 13 min. Optimum polymerization takes place in the pressure pot at 45 °C at 2 bar (10 min). If processing without pressure curing unit, apply large amounts of MEGA-A in layers. Each layer must be completely cured and chilled before applying new dough.

### *-Preparation and polishing*

After polymerization, the denture is removed from the model and processed with cross-cut milling cutters. Before final polishing with pumice stone, the denture is smoothed with abrasive paper of decreasing grain size. After production, MEGA-A contains an admissible residual amount of monomer < 2 % by weight. In order to prevent the risk of incompatibilities further, it is recommended to store the prosthesis in lukewarm water for at least 12 hours before it is inserted.

### *- Residual monomer content*

The residual amount of monomer after polymerisation is <4,5 % by weight. In order to prevent the risk of incompatibilities further, it is recommended to store the prosthesis in lukewarm water for at least 12 hours before it is inserted.

## **9. Shelf life**

Liquid: 2 years

Powder: 3 years

## **10. Precaution and storage**

Not above 25 °C

## **11. Delivery units**

Powder: 80 g / 500 g / 1000 g

Liquid: 80 ml / 500 ml / 1000 ml

**12. Waste treatment**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation. Handle contaminated packages in the same way as the substance itself.

**13. Date information prepared**

2018-08-30

**14. Manufacturer and address**

MEGADENTA Dentalprodukte GmbH  
Carl-Eschebach-Str. 1A / D-01454 Radeberg, Germany

---

1200 0818 531