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Shades: A 1; A2; A3; A3,5; A4 ; B1 ; B2 ; Dentin A2 ; Dentin A3

2. Composition

The inorganic filler particles comprise barium aluminium bore glass and highly dispersed silicon dioxide. The monomer matrix contains:

2,2-Bis-4(2hydroxy-3-methacryloxy-propyloxy)-phenyl-propane,

7,9,9-trimethyl-4,13-dioxo-3,14-dioxa-5,12-diaza-hexadecan-1,16-di-yl-bis-methacrylate,

1,4-Butandioldimethacrylate.

3. Indication

Usage in dentistry for Class V restorations (cervical caries, root erosion, wedge-shaped defects), small posterior restorations, anterior restorations (classes III and VI) restorative therapy for mini-cavities, preventive resin restorations in premolars and molars, adhesive cementation of ceramic and composite restorations.

4. Contraindication

- If a dry working field cannot be established.
- If the patient is known to be allergic to any of the ingredients.

5. Side effects

The unpolymerized material may cause contact allergy.

6. Other interactions

Megafill Flow may not be used together with underfillings containing eugenol because phenolic substances inhibit the polymerisation process.

7. Warning

Megafill Flow is light sensitive as polymerization can happen with day-light. Therefore the light exposure should be limited (take material from syringe just before application, close syringe directly after that).

8. Hazard Information

Read Safety data sheet first!

General information

Polymerized material is non-hazardous.

After inhalation

Provide fresh air. When in doubt or if symptoms are observed, get medical advice.

After contact with skin

After contact with skin, wash immediately with polyethylene glycol, followed by plenty of water. Take off immediately all contaminated clothing and wash it before reuse. Medical treatment necessary.

After contact with eyes

In case of contact with eyes flush immediately with plenty of flowing water for 10 to 15 minutes holding eyelids apart and consult an ophthalmologist.

After ingestion

Rinse mouth immediately and drink plenty of water.

Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

9. Application

-Cleaning, shade selection and drying the teeth:

- Before etching the tooth surface has to be thoroughly cleaned with an oil free polishing paste or a pumice powder.
- Shade matching should preferably be carried out in daylight using the Vita classic shades
- Dry the surfaces.
- Use of rubber dam is recommended.

-Enamel etching:

- Use the supplied brush to etch the cleaned enamel with C-Cid.
- Rinse C-Cid with water after 20-40s thoroughly and dry carefully with oil- and water-free air. The etched and dried enamel must show a smooth, chalky surface.
- In case of additional Dentin etching (total-etch technique), allows C-Cid to act on the Dentin surface for max. 15 s, rinse off with water and slightly fade with air
- The etched enamel must not be contaminated by saliva or blood since this would lead to a deposit of protein, which would inactivate the etched enamel surface. If this happens, repeat the etching process.

-Application:

- The material is applied with the needle directly into the cavity.
- Remove any air bubbles with the needle or with a suitable instrument.
- Polymerize with a light source (e.g. Megalux) for 40 s.
- Restorations deeper than 2,5 mm should be filled in layers, each layer have to be cured separately (multi-phase technique).
- The final forming and the removal of excess composite are carried out with rotating finishing instruments. For contouring and finishing fine and ultra fine finishing diamonds are used. For polishing flexible discs and polishing strips are suitable

10. Shelf life

3 years

11. Precaution and Storage

We recommend storage in a fridge (+5 °C [41 °F] to +8 °C [46 °F]); do not exceed 25 °C [77 °F].

12. Delivery units

single syringe Megfill Flow á 2 g

assortment 4 x Megafill Flow á 2 g

Accessories: brush holder, brush inserts, Luer-Lock cannulas

13. Waste treatment

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

Keep out of the reach of children!
For use in dentistry only!

14. Date information prepared

2018-10-18, Rev. 3

15. Manufacturer and Sales

MEGADENTA Dentalprodukte GmbH

Carl-Eschebach-Strasse 1 A

D-01454 Radeberg, Germany