Megaseal N

Basis-UDI-DI: + EMD12VLACKNMX

[de] Gebrauchsanweisung
[en] Instructions for use

 0124

MEGADENTA
Dentalprodukte

MEGADENTA Dentalprodukte GmbH
Carl-Eschebach-Str. 1A, D-01454 Radeberg, Germany
Telefon +49(0)3528 453-0, Fax +49(0)3528 453-21
E-Mail: info@megadenta.de www.megadenta.de

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[en] Instructions for use

Read and keep before using the product! Medical device for professionals only

1 General safety instructions

Observe safety data sheet!

P210 Keep away from heat. No Smoking.

H225 Highly flammable liquid and vapor.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

H335 May cause respiratory irritation.

H360Fd May damage fertility. Suspected of damaging the unborn child.

H412 Harmful to aquatic life with long lasting effects.

Caution! Megaseal N contains methyl methacrylate! Methyl methacrylates and their vapor are highly flammable. Do not smoke, keep away from sources of ignition!

Boiling point: 100 °C, Ignition temperature: 430 °C

Methyl methacrylates cause skin irritation, can cause allergic skin reactions and irritate the respiratory tract.

Warning! Megaseal N contains the photoinitiator TPO (diphenyl(2,4,6-trimethylbenzoyl) phosphine oxide), a hazardous substance classified as a CMR substance according to Regulation (EC) No. 1272/2008. When handling, the protective measures listed in the safety data sheet and below must be strictly observed. When light polymerisation is carried out correctly, TPO is largely converted; residual levels in the polymerised product are very low and, according to current knowledge, are within a toxicologically safe range..

Avoid body contact with unpolymerized product! Never use unpolymerized product directly on the patient (only extraorally / outside the oral cavity)!

Use protective equipment: Protective goggles/face protection, suitable protective gloves with protective function, especially against methacrylates (e.g. polyethylene gloves, 4h gloves); change contaminated clothing.

After skin contact: Rinse immediately with plenty of water. Medical treatment necessary if symptoms occur (irritation of the skin).

After eye contact: Rinse immediately with running water for 10 to 15 minutes with the eyelids open and consult an ophthalmologist.

Wear respiratory protection in case of insufficient ventilation.

The nanomaterial contained in the polymerized state does not cause any increased risk regarded

to cell-permeable nanoparticles, since the polymerization results in a solid compound.

During the period of use of a sealed prosthetic material, which may be several years, the leaching of chemicals/substances from the polymer matrix cannot be completely excluded. This depends on various factors (e.g. the patient's dietary habits). The patient must be informed about the residual risk of genotoxic, cytotoxic and/or sub-acute systemic toxic effects, with a particular focus on pregnant women, nursing mothers and children.

All serious incidents related to the product must be reported to the manufacturer and the competent authority of the member state where the user and/or patient is established. Please direct any feedback regarding the product to: reclamation@megadenta.de

2 Description

2.1. General description

Megaseal N is a transparent, light-curing, single-component varnish for sealing dental prosthetic surfaces such as denture bases, temporary crowns and bridges, or for impression trays made of light-curing tray material. Megaseal N has a slightly watery consistency and is therefore easy to spread. It seals and smooths the surfaces and thus protects against discoloration and delays the accumulation of plaque.

2.2. Composition

Megaseal consists of a mixture of acrylates and methacrylates with integrated nanoparticles and contains additives such as surfactants, stabilizers, photoinitiators.

3 Intended use

3.1. Warnings

The liquid and its vapor are highly flammable. Remove all sources of ignition. Close the bottle immediately after each use.

Megaseal is sensitive to light, so that polymerization can be triggered by daylight. Removal from the primary packaging (bottle) should therefore only be made directly before use and the packaging should be resealed immediately.

If the primary packaging is damaged, the product must be disposed of in accordance with 3.12.

3.2. Intended purpose / Intended use

Megaseal is intended for sealing dental prosthetic acrylic materials.

Megaseal can also be used for sealing other similar acrylic surfaces such as impression trays made of light-curing tray material.

The user is obliged to use the material in accordance with the intended purpose specified by the manufacturer MEGADENTA Dentalprodukte GmbH. The manufacturer accepts no liability for damage resulting from use for other purposes.

3.3. Clinical benefit

The clinical benefit of Megaseal lies in the smoothing of the surfaces of acrylic prostheses, thereby reducing plaque accumulation, discoloration, and mechanical irritation of the oral mucosa while significantly improving the aesthetics and functionality of the dentures.

3.4. Indications

Smoothing of dental prosthetic acrylic surfaces

- to reduce mechanical irritation of the tongue and mucous membrane
- to reduce plaque accumulation
- to reduce discoloration
- to maintain an aesthetic appearance.

3.5. Contraindication

Known hypersensitivity to the substances listed in the composition, especially methacrylates

3.6. User profile

Megaseal is intended for professional use by dental technicians. It is intended for use in the dental laboratory or dental practice laboratory. No special training is required.

3.7. Patient group

The use of Megaseal is not restricted to any particular patient group. Contraindications and safety instructions must be observed.

3.8. Undesirable side effects

Megaseal contains methacrylates. Contact allergies may occur in individual cases.

3.9. Undesirable interactions

No undesirable interactions are known.

3.10. Special storage instructions

Keep bottle tightly closed. Store away from light, not above 25 °C and in a well-ventilated place. Storage outside the specified conditions may result in damage to the product, e.g. deterioration in mechanical properties due to changes in composition caused by evaporation of components, premature curing. It can then no longer be used as intended and must be disposed of in accordance with 3.12.

3.11. Shelf life

3 years.

Do not use after the expiry date stated on the label!

3.12. Waste treatment

This product and its container must be disposed of as hazardous waste. Contaminated packaging must be treated in the same way as the substance. Do not allow product or primary packaging that is not completely empty to enter drains, waterways, subsoil or soil. Dispose of in accordance with official regulations.

Non-contaminated and completely emptied packaging as well as fully cured product residues can be disposed of with normal municipal waste (household waste). Observe waste separation.

4 Application

The surface of the denture to be sealed should be as smooth as possible and free of irregularities. It must be cleaned and dried before sealing. The surface must be free of oils and grease.

Megaseal N is applied evenly to the plastic surface with a brush in a thin film. After a short drying time (evaporation time of the solvent approx. 30 s), the wetted surface is polymerized in a commercially available light polymerization device (e.g. Megalight) with UV or blue light in the light wave range between 320 - 500 nm for approx. 3-5 min. The polymerization time depends on the intensity of the device. High-power devices (e.g. stroboscopic devices) enable a short polymerization time. However, this should not be less than 3 minutes. The instructions for use of the device used must be observed! The prosthesis must be positioned in the device so that the surfaces to be sealed are completely irradiated; if necessary, rotate and/or turn and repeat the irradiation.


MMA/PMMA-based materials always contain a small amount of remaining monomer after manufacture. To reduce the risk of intolerance reactions, the coated prosthesis, crown or bridge should be stored in lukewarm water for at least 12 hours before insertion. No other finishing work is necessary.

5 Delivery units

Article no. 2VLACK-N 20 ml bottle with brush insert
Article no. 2VLACK-N3 30 ml bottle with drip insert

6 Symbols used in product labeling

	Manufacturer
	Medical device
	Item number
	Batch

	Usable until
	Protect from sunlight
	Upper temperature limit
	Observe instructions for use
	Irritant
	Highly flammable liquid and vapor
	Hazardous
	Harmful to aquatic life
	Contains nano materials
	Contains hazardous substances (see chapter 1)
Rx only	Prescription drug (US Federal Law restricts this to sale by or on the order of a dentist)
	European Conformity Mark

1690 0825 538 / 2025-08-21 (Rev. 4)