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Instruction for use (GB)

1. Product decription

Megaseal Universal is suitable for use as a light-curing varnish for sealing of acrylic dentures (PMMA), temporary crowns and bridges as well as tray plates made of light-curing tray material. The material is intended exclusively for use in the dental laboratory or for extraoral application (outside the oral cavity) in the dental practice.
Megaseal Universal seals the surfaces and thus serves to protect against discoloration and delays the accumulation of plaque.

2. Application

Megaseal Universal is applied to the clean & dry plastic surface as a thin film using a brush. Megaseal Universal has a viscosity like water and is therefore easy to spread. After a short drying time (solvent evaporation time approx. 30 s), the wetted surface is polymerized in a commercially available light polymerization unit with UV- or blue light in the light wave range between 320–500 nm (e.g. Megalight unit) for approx. 3–5 min. The polymerization time depends on the intensity of the device. High-performance devices (e.g. stroboscopes) allow a short polymerization time. However, this should not be less than 3 minutes.

3. Composition / Technical datas

The light-curing surface sealer consists of a mixture of acrylates and methacrylates.
Composition: Multifunctional methyl methacrylates, acrylates, photoinitiators, stabilizers.

4. Importend advice

Attention! Megaseal Universal contains methyl methacrylates and is flammable!
Boiling Point: 100°C
Ignition temperature: 430°C
Methyl methacrylate and its vapor are highly flammable, cause skin irritation and can cause allergic skin reactions. It can irritate the respiratory tract.
Do not smoke, keep away from sources of ignition. Avoid skin contact in the unpolymerized state. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
Wear goggles/face protection when working. In case of contact with skin, wash immediately with plenty of water.
Wear respiratory protection if there is insufficient ventilation.

5. Shelf life

3 years

6. Storage and waste conditions

Do not store above 25°C! Megaseal Universal contains methyl methacrylates and is flammable!
Store in a well-ventilated place. This product and its container must be disposed of as hazardous waste. Contaminated packaging must be treated in the same way as the substance.

7. Package and delivery units

20 ml bottle with brush insert	REF: 2VLACK-UV2
30 ml bottle with drip insert	REF: 2VLACK-UV

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