# Megatray®

[de]	Gebrauchsanweisung
[en]	Instructions for use
[bg]	Инструкции за употреба
[cs]	Návod k použití
[da]	Brugsanvisning
[es]	Instrucciones de uso
[fi]	Käyttöohje
[fr]	Mode d'emploi
[hr]	Upute za uporabu
[hu]	Használati utasítás
[it]	Istruzioni per l'uso
[lt]	Naudojimo instrukcija
[lv]	Lietošanas pamācība
[nl]	Gebruiksaanwijzing
[no]	Bruksanvisning
[pt]	Instruções de utilização
[ro]	Instrucțiuni de utilizare
[ru]	Инструкция по применению
[sk]	Návod na použitie
[sv]	Bruksanvisning

**Basis-UDI-DI:**  
**+EMD11MYDVOKS2 (Megatray)**  
**+EMD11MYDVOKRS125 (Megatray Basisplatte)**

**MEGADENTA**

Dentalprodukte

MEGADENTA Dentalprodukte GmbH  
Carl-Eschebach-Str. 1A, D-01454 Radeberg, Germany  
Telefon +49(0)3528 453-0, Fax +49(0)3528 453-21  
E-Mail: [info@megadenta.de](mailto:info@megadenta.de) / [www.megadenta.de](http://www.megadenta.de)

Produktname / Product name	Varianten / Variants (colors)	Einheit / Unit	Artikel-Nr. / Item no.
Megatray	Einzelpack blau / Single pack blue	1 Dose á 50 Stück / 1 box of 50 pieces	1MYDVOK-BL
	Einzelpack weiß / Single pack white		1MYDVOK-WS
	Einzelpack transparent / Single pack transparent		1MYDVOK- TR
	Einzelpack rosa / Single pack pink		1MYDVOK-RS
	Doppelpack blau / Double pack blue	2 Dosen á 50 Stück / 2 boxes of 50 pieces	1MYEVOK-BL
	Doppelpack weiß / Double pack white		1MYEVOK- WS
	Doppelpack transparent / Double pack transparent		1MYEVOK-TR
	Doppelpack rosa / Double pack pink		1MYEVOK-RS
Megatray Basisplatte / Megatray Base plate	Einzelpack rosa / Single pack pink	1 Dose á 50 Stück / 1 box of 50 pieces	1MYDVOK-RS1
	Doppelpack rosa / Double pack pink	2 Dosen á 50 Stück / 2 boxes of 50 pieces	1MYEVOK-RS1
	Start-Kit rosa / Start-Kit pink	1 Dose á 25 Stück / 1 box with 25 pieces	1MYDVOK-RS2

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## [en] Instructions for use

Read and keep before using the product! Medical device for professionals only

### 1 General safety instructions

Observe safety data sheet!

H317: May cause an allergic skin reaction. Avoid contact of uncured material with skin/mucous membrane and eyes. It may be slightly irritating and in rare cases may cause sensitization to methacrylates.

Use protective equipment: Safety glasses, suitable protective gloves with protection function to Methacrylates (e.g. Polyethylene-gloves, 4h-gloves). Do not wear contaminated work clothing outside the workplace.

After skin contact: Wash immediately with plenty of soap and water. Medical treatment necessary if symptoms occur (irritation of the skin).

After eye contact: Immediately flush with running water for 10 to 15 minutes with the eyelid open and consult an ophthalmologist.

When grinding the polymerized sheets, use a suitable dust extractor and follow the instructions for use of the abrasive manufacturer, especially with regard to dust protection. If necessary, use a self-contained breathing apparatus.

All serious incidents related to the product must be reported to the manufacturer and the competent authority of the member state where the user and/or patient is established. Please direct any feedback regarding the product to: [reclamation@megadenta.de](mailto:reclamation@megadenta.de)

## 2 Description

### 2.1 General description

Megatray is a pre-shaped, modelable, light-curing composite for the production of individual impression trays in the dental laboratory or dental practice. It features low shrinkage, low water absorption and high dimensional stability.

Available variants:

- Megatray: plate thickness approx. 2.2 mm in the colors transparent, blue, pink and white
- Megatray base plates: plate thickness approx. 1.4 mm in pink

### 2.2 Composition

Acrylates, methacrylates, HQME, fumed silica, glass beads, photochemical initiators, pigment, vaseline (packaging aid).

## 3 Intended use

### 3.1 Warnings

Megatray is light-sensitive, so that polymerization can already be triggered by exposure to daylight. Therefore, the material should not be exposed to light unnecessarily (removal from the box only immediately before adaptation to the model, immediate sealing of the box after removal). In case of damage to the primary packaging (box), dispose of the product according to 3.12.

### 3.2 Intended purpose / intended use

Megatray is intended for dental use:

- for the production of individual impression trays for precision impressions for accurately fitting restorations with the highest marginal accuracy, e.g.
  - functional impressions
  - situational impressions for unfavorable dental arch and jaw shapes and/or deep-set bands

Megatray base plates are intended for dental use:

- for production of a occlusal base plate / for intraoral support pin registration

The user is obliged to use the material in accordance with the intended purpose specified by the manufacturer MEGADENTA Dentalprodukte GmbH. The manufacturer accepts no liability for damage resulting from use for other purposes.

### 3.3 Clinical benefit

The clinical benefit of light-curing impression materials lies indirectly in the replacement of the dental apparatus and thus in a restoration of masticatory function, since precision and/or functional impressions are essential for the accurately fitting production of dentures.

### 3.4 Indications

Megatray:

- Precision impressions for accurately fitting restorations with the highest marginal accuracy, such as
  - functional impressions
  - Situation impressions for unfavorable dental arch and jaw shapes and/or deep-set bands

Megatray base plates:

- Bite registrations / Support pin registrations

### 3.5 Contraindication

In case of known hypersensitivity to methacrylates, skin contact should be avoided.

### 3.6 User profile

Megatray is a product for professional use by dental technicians or dentist. It is intended for use in the dental laboratory or dental practice. No special training is required.

The finished products made of the material (individual impression trays, base plates, support pin registrations) are intended for professional use by dentists in a dental practice. Special training is not required.

### 3.7 Patient group

The use of Megatray is not limited to any particular patient group. Contraindications must be observed.

### 3.8 Undesirable side effects

Megatray contains methacrylates. These may cause allergic skin reactions. In rare cases sensitization is possible.

### 3.9 Undesirable interactions

No interactions are known.

### 3.10 Special storage instructions

Keep the box tightly closed. Store protected from light and at 4-25°C.

Storage outside the specified conditions may cause damage to the product, e.g. premature curing, stickiness. It can then no longer be used as intended and must be disposed of in accordance with 3.12.

### 3.11 Shelf life

3 years from manufacture. Do not use after the expiry date stated on the label!

### 3.12 Waste treatment

Do not allow product or primary packaging which has not been emptied of all residues to reach sewage system, water, subsoil or soil. Disposal according to official regulations.

Non-contaminated and completely empty packaging as well as fully cured products can be disposed of with normal municipal waste (household waste). Observe waste separation.

## 4 Application

- The situation model should be prepared as required and insulated with an alginate insulation (e.g. Mega-I). Undercutting areas of the model should be blocked out before.
- Open the box at the handle tab, remove one Megatray plate incl. protective fleece and immediately reseal the box tightly.
- The protective fleece is removed from the Megatray sheet.
- The sheet is adapted to the insulated model and trimmed at the edges.
- The handle can be formed from the excess material. The handle is not used when making a base plate.
- The sheet is then polymerized on the model in a blue light or UVA light device (e.g. Megalight) for 5 minutes. The instructions for use of the light unit used must be observed!
- The sheet is removed from the model and polymerized from the back for another 5 minutes.
- The individual impression tray made in this way can be processed immediately.
- The inhibition layer (adhesive layer on the top of the material) allows subsequent correction (e.g. reinforcement of the sheet). The inhibition layer can be easily removed with a moistened tissue with an alcoholic cleaning agent (e.g. Megaclean).



## 5 Delivery units

Megatray: 1 box of 50 pieces / 2 boxes of 50 pieces

Megatray Base plate: 1 box of 50 pieces / 2 boxes of 50 pieces / 1 box with 25 pieces

## 6 Symbols used in product labeling

	Manufacturer
	Medical device
	Item number
	Batch
	Usable until
	Protect from sunlight

	Temperature limit
	May cause an allergic skin reaction
	Observe instructions for use
	Prescription drug (US Federal Law restricts this to sale by or on the order of a dentist)
	European Conformity Mark

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