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Directions of use

N-Fill Flow

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Directions of use

1. Information

N-Fill Flow , light-curing flowable nano concept filling material
Shades: A1; A2; A3; A3,5; B2

2. Composition

The inorganic filler particles comprise strontium aluminium bore glass, high-dispersal non-agglomerated silicon dioxide and agglomerated silicon dioxide.

The monomer matrix contains:

2,2-Bis-4(2hydroxy-3-methacryloxy-propyloxy)-phenyl-propane,

1,4-Butandioldimethacrylate,

Triethyleneglycol-dimethacrylate.

3. Indications

Class V restorations, small restorations class I and II, anterior fillings of class III and IV
(cervical caries, root erosion, wedge-shaped defects, posterior restorations, anterior restorations)

4. Side Effects

Placement of N-Fill Flow is contraindicated if the patient is known to be allergic to any of the ingredients in N-Fill Flow .

5. Interaction

Application of eugenol containing material, specifically base liners, is not indicated, it will impair the polymerization.

6. Warning

The C-Cid etching contains 35% phosphoric acid. Avoid skin- and eyes contact with C-Cid and unpolymerized material. N-Fill Flow cures with light, it is sensitive to blue light (operating light, sun light). For that reason bright light during application of N-Fill Flow should be avoided.

7. Application

-Shade Selection

The tooth and the adjacent distal surfaces has to be cleaned with a bristle brush and prophylactic paste. Shade matching should preferably be carried out in daylight using the MEGADENTA Shade Guide. The shades corresponding to the Vita Shade[®] system.

-Cavity Preparation

If the working field cannot be established dry the application will be contraindicated. The cavity is prepared and dried in accordance with the general rules of composite filling treatment, unprepared dentine is cleaned with a cleaning paste. Caries-free cervical defects are not prepared, but merely cleaned with pumice and a suitable cleaning paste using a rubber cup or rotary brush. The prepared cavity is now rinsed with water and dried with oil free air.

-Base Liner

To protect exposed dentine the application of a glass ionomer liner, phosphate cement or calcium hydroxide base liner (deep cavities) is recommended.

-Acid Etching

The cleaned enamel edges are etched with C-Cid for 20 – 40 s. After 20 - 30 s C-Cid is applied additional on the dentine allowed to take effect for another 10 - 15 s and then thoroughly rinsed with water. Excessive drying should be avoided in this case, as this would cause the exposed collagen fibre to collapse.

-Application of the Primer

The Primer has now to be used after the recommendations of the producer. Do not rinse off the dentin primer! The cavity is dried thoroughly with oil and water free air or until the dentine surface appears dull and dry. It is essential that the primed dentine and the etched enamel are dry and contaminant free for the bond application. C-Bond is applied onto the primed dentine in the dried cavity and the etched enamel, blown out thinly with oil free air and then cured for 10 s with a blue light source (e.g. Megalux).

-Application of N-Fill Flow

After the polymerization of C-Bond N-Fill Flow is applied. A direct application is possible with the syringe needles enclosure. To avoid an oxygen inhibited layer, a coat of oxygen-tight glycerine gel can now be applied along the restoration margins. Make sure that the glycerine gel is not mixed with the N-Fill Flow.

After completed contouring the restoration is cured for 20-40 s with a blue light source (e.a. Megalux). The beam has to be directed to all surfaces of the restoration. Restorations deeper than 2,0 mm should be filled in layers, each layer have to be cured separately (multi-phase technique).

-Shaping and Finishing

The final forming and the removal of excess composite are carried out with rotating finishing instruments. For contouring and finishing fine and ultra fine finishing diamonds are used. For polishing flexible discs and polishing strips are suitable. All shaping and finishing have to be carried out under tooth preserving and with water-cooling.

8. Shelf life

3 Years

9. Precautions and Storage

Not to be stored above 25°C

Storage in refrigerator (+5 to +8°C) prolongs shelf life

10. Deliver units

syringe of 2,0g

11. Date information prepared

January 2007

12. Manufacture and Sales

MEGADENTA Dentalprodukte GmbH

D-01454 Radeberg, Germany